




**Title: Technical Datasheet for Disposable Medical Face Mask – SM120**

Product Name	Disposable Non-Sterile Medical Face Masks Type IIR 3 Ply
Product Code	SM120
Product Branding	MAX
Product Size	180mm x 90mm (ear loops= 180mm)
Product Weight	2g
Product Material	Inner layer: White spunbonded non-woven polypropylene 25g/sqm Outer layer: Blue spunbonded non-woven polypropylene 25g/sqm Filter media: Meltblown white non-woven polypropylene 25g/sqm, 0.35mm thickness Elastic strap: 2mm thick white elastic – latex free Nose piece: 3mm thick metal wire encased in white plastic, 165mm length
Product Performance	<u>EN14683:2019</u> Bacterial Filtration Efficiency (BFE) >= 99% Differential Pressure <= 10.2 Pa/cm <sup>2</sup> Synthetic Blood Resistance: PASS Microbial Contamination/Bioburden <= 23 C.F.U./g
Standard Product Colour	Blue/White
Packaging Size	70mm x 200mm x 100mm
Packaging Units	50 per inner box – 40 boxes per outer box
Manufacturer	CMT Equipment LTD – Trident Works, Mulberry Way, Belvedere, London, DA17 6AN
Storage Conditions	Store in a well-ventilated room where the relative humidity does not exceed 80%. Avoid high temperatures, keep away from fire and flammable materials.
Shelf Life	5 years in original packaging
Handling	Keep in original packaging until use – do not reuse. Replace when damp or damaged
Conformities	EN 14683:2019; CE Marked; Registered with MHRA. Produced in accordance with MDR 2017/745. Class I medical device (non-sterile)
Disposal Instructions	Dispose of contents in accordance with local regulations. For used masks, follow medical advice and regulations for disposal of biohazardous items if relevant.



**Title: Technical Datasheet for Disposable Medical Face Mask – SM120**

Usage Instructions	<p style="text-align: center;"><b>INSTRUCTIONS</b></p> <div style="display: flex; justify-content: space-around; text-align: center;"> <div style="border: 1px solid black; padding: 5px; width: 15%;"> <p><b>1</b></p> <p>Hold the mask to your face to covering your nose bridge</p> </div> <div style="border: 1px solid black; padding: 5px; width: 15%;"> <p><b>2</b></p> <p>Stretch each strap over your ears so it's comfortable</p> </div> <div style="border: 1px solid black; padding: 5px; width: 15%;"> <p><b>3</b></p> <p>Adjust the upper and lower parts to fit your face</p> </div> <div style="border: 1px solid black; padding: 5px; width: 15%;"> <p><b>4</b></p> <p>Use both fingers to press in the metallic strip</p> </div> </div> <p style="text-align: center;"><b>PRECAUTIONS:</b></p> <p>One for one-time use. Do not use if package is damaged, destroy if so.          Use with caution if you react to any nonwoven fabrics. Does not contain latex.          Store in a well ventilated room where the relative humidity does not exceed 80%.          Avoid high temperatures, keep away from fire and flammable materials. Dispose of contents in accordance with local regulations. For used masks, follow medical advice and regulations for disposal of biohazardous items if relevant. Do not reuse.</p>
Dermatological Tolerance	No known intolerance, does not contain latex
Conditions	Non-sterile
Product Image	

Packaging Image

